PRJ 22

R1
Cytotoxic activity of Bio-normalizer against human cancer cell lines
de Castro-Bernas G, Remo-Husmillo G, Osato JA,
Santiago LA and Mori A

Cytotoxic activity of Bio-Normalizer against human cancer cell lines

Gloria de Castro-Bernas¹, Gemma Remo-Husmillo², James Akira Osato^{3,4}, Librado A. Santiago^{3,5} and Akitane Mori⁵

¹Research Center for the Natural Sciences, University of Santo Tomas, España, Manila, Philippines, ²Osato Research Foundation, Inc., 12th Floor, Royal Match Building, 6780 Ayala Avenue, Makati Metro Manila, Philippines, ³Osato Research Institute, 4th Floor Daiichi Taiyo Bldg. 5-567-9 Tonomachi, Gifu City, Gifu 500, Japan, ⁴The United Graduate School of Agricultural Science, Gifu University, 1-1 Yanagido, Gifu 500, Japan, ⁵Department of Neuroscience, Institute of Molecular and Cellular Medicine, Okayama University Medical School, 2-5-1 Shikata-cho, Okayama 700, Japan

Address correspondence to:

Librado A. Santiago, D.M.Sc. Osato Research Institute 4th Floor Daiichi Taiyo Bldg. 5 - 567 - 9 Tonomachi, Gifu City Gifu 500, Japan Telefax: 81-(58)-247-4496

Abstract

The cytotoxicity of Bio-Normalizer, a fermented nutritional supplement, against four human cancer cell lines was determined using the tetrazolium-based colorimetric assay. Based on the IC₅₀ values derived, epidermoid larynx carcinoma (HEp-2) was the most sensitive to Bio-Normalizer. Bronchio-alveolar adenocarcinoma (SLNI-52) and lung adenocarcinoma (A 549) showed comparable sensitivity to Bio-Normalizer while a higher IC₅₀ value was obtained for cervical epitheloid carcinoma (HeLa).

A separate study demonstrated that at a reduced concentration Bio-Normalizer can increase cell viability of normal human foreskin fibroblasts to more than 100%.

These findings suggest a potential cytotoxic activity for Bio-Normalizer against the carcinoma cells used.

Keywords: cytotoxicity, Bio-Normalizer, cancer cells

1. Introduction

To scientifically document Bio-Normalizer's purported therapeutic effect on cancer cases its potential cytotoxic activity against human cancer cell lines was determined. Human carcinoma cell lines are useful models for in vitro assessment of potential cell killing effects of commercial drugs or natural products. Epithelial cells are of particular interest because most adult human cancers are carcinomas [1].

Bio-catalyzer α . ρ No. 11 (Bio-Normalizer), a commercial natural fermented health food supplement marketed in Japan and in the Philippines has proven scavenging action on free radicals [2] and was found to increase superoxide dismutase activity in the brain of aged rodents as a defense mechanism against age-related accumulation of superoxide radicals [3]. Recently, it was observed to enhance gamma interferon in the body [4] advocating its significant role in the obstruction of diseases. Its antigenotoxic potential against known mutagens and carcinogens like dimethylnitrosoamine and cyclophosphamide (in vivo) has been established [5].

The use of Bio-Normalizer was observed to have improved various medical cases usually mediated by free radicals. Free radicals damage biomolecules like proteins and DNA, causing diseases such as cancer [6], cardiovascular diseases [7] and brain dysfunctions [8]. Radiation-induced carcinogenesis may be initiated by free radicals [9]. Free radical intermediates, in turn, damage DNA. Base hydroxylation and strand breaks occur leading to a fall in ATP levels and death of the cell. In surviving cells, mutation takes place leading to malignant transformation as evident in naturally occurring human cancers of the colon, pancreas and lung [10,11]. Cancer is characterized by an increase in the number of abnormal cells, their invasion of adjacent tissues, and lymphatic spread of malignant cells to distant sites [12].

2. Materials and Methods

2.1 Preparation of compounds

Bio-Normalizer (Sun-O International, Inc., Gifu, Japan) was dissolved in phosphate buffered saline (PBS) with concentrations ranging from 20 - 100 mg/ml while the positive control, vincristine sulfate (Pharmachemie B.V., Haarlem, Holland), was obtained as formulated for clinical use and reconstituted with PBS to attain a 100 μg/ml stock solution. They were prepared freshly for each experiment and finally passed through a 0.22 μm sterile filter prior to use. Trypsin and RPMI 1640 supplemented with 0.03% L-glutamine and 0.20% sodium bicarbonate were purchased from Sigma Chemical Co. (St. Louis, MO. USA).

2.2 Cell cultures

Human epidermoid carcinoma of the larynx cell line (HEp-2) and epitheloid cervix carcinoma cell line (HeLa) were obtained from the Research Institute for Tropical Medicine (Alabang, Metro Manila, Philippines). Human lung adenocarcinoma cell line (A 549) and bronchio-alveolar adenocarcinoma cell line (SLNI-52) were gifts from the Section of Hematology and Oncology, University of California at Los Angeles (Sepulveda, California, USA). Cells were cultured in RPMI 1640 supplemented with 0.03% L-glutamine, 0.20% sodium bicarbonate, 10% calf bovine serum (GIBCO, New York, USA), penicillin, streptomycin and fungizone (GIBCO, New York, USA) at 37°C in a humidified atmosphere with 5% CO2.

2.3 MTT assay

A modification of the quantitative colorimetric assay for mammalian cell survival and cell proliferation by Mossman [13] using 3-(4,5-dimethylthazol-2-yl)-2,5-diphenyl tetrazolium bromide (MTT) was adapted in this study.

In the preliminary experiments done, the optimal seeding density required for each cell line to be in the log phase when the assay was terminated 4 days after plating was determined to be 10,000 cells per well. To the overnight cultures of cells, Bio-Normalizer (20 - 100 mg/ml), vincristine sulfate (20 μ g/ml) and equivalent amount of PBS were added. Each

concentration point for the test compounds was representative of 4 - 8 multiplicate wells for each cell line. The cells were incubated at 37°C for 4 days and then 125 μ g MTT (5 mg/ml stock solution) was added to each well. For the solubilization of the formazan crystals, modification by M.B. Hansen, et al [14] was followed. One hundred microliter of solubilizing solution (20% sodium dodecyl sulfate in 50% dimethylformamide, pH 7.4) was added to each well and incubated overnight at 37°C. Optical density of each well was measured manually at 570 nm using Photometer 4010 (Boehringer Mannheim, Germany).

Response curves were plotted, with the IC_{50} values for each cell line determined graphically as the dose of Bio-Normalizer causing a 50% reduction in absorbance as compared to the control values.

2.5 Statistical analysis

The results obtained were analyzed by the student's t-test using the statistical significance at the 0.05 probability level.

3. Results

Tumor cell killing as a function of Bio-Normalizer concentration was determined for the four available cell lines and this was repeated in 2-3 experiments. Bio-Normalizer was tested over a concentration range of 20-100 mg/ml.

Results showed significant differences in the effects of Bio-Normalizer on the four cell lines (p < 0.05). The IC₅₀ value of Bio-Normalizer for HEp-2 was calculated to be 40 mg/ml; 60 mg/ml for A 549 and 64 mg/ml for SLNI-52. The IC₅₀ for HeLa was found to be the highest at 92 mg/ml. The cell killing effect of Bio-Normalizer on A 549 and SLNI-52 was observed to be comparable to the activity of vincristine sulfate at 20 μ g/ml (p > 0.05) while the solvent used to dissolve the test compounds did not have cytotoxic effects on the four cell lines used (Figure 1).

It was observed that Bio-Normalizer kills cancer cells in a dose-dependent manner. At the highest concentration of Bio-Normalizer (100 mg/ml), 22, 26, 30, and 50 per cent cell viability was observed on HEp-2, SLNI-52, A 549, and HeLa, respectively, while at its lowest concentration (20 mg/ml) cell viability was 90, 78, 82, and 84 per cent, respectively (Figure 2). Bio-Normalizer starts killing cancer cells at about 20 mg/ml.

Microscopic observations showed that long exposure of the cancer cells to Bio-Normalizer (4 days prior to assay) lead to intracellular granulations, disintegration of the cell membrane, and inability of the cells to adhere to the substratum and proliferate (data not shown).

In a separate study, it was found out that at 10 - 55 mg/ml and 100 mg/ml Bio-Normalizer, cell viability of the normal human foreskin fibroblasts (HFF) was more than 100% and 70%, respectively. It starts killing HFF at 58 mg/ml.

4. Discussion

The heterogeneity of response of Bio-Normalizer on different tumor cell populations suggests not only a cytotoxic activity but anti-neoplastic activity as well. Observations on the action of Bio-Normalizer were made on actively dividing cell population [15] which distinguishes an anti-neoplastic activity from a non-specific cytotoxic activity.

The uncontrolled growth of some neoplastic cells can be attributed to the concept of autocrine production of growth factors [1]. Ectopic hormones like bombesin and human chorionic gonadotropin produced by carcinomas are examples of autocrine growth factors. These hormones enhance the growth of epithelial cells *in vitro* by binding to specific membrane receptors and elicit cascades of reactions leading to cell growth and differentiation. Bio-Normalizer could have inhibited the affinity of these hormones to the membrane receptors, thereby, inducing irreversible inhibition of DNA and membrane protein synthesis.

In a previous *in vivo* study, it has been shown that Bio-Normalizer inhibits the toxohormone-L-induced lypolysis [16] and increases significantly the natural killer cell activities of tumor bearing mice [17]. Natural killer cells were able to recognize and lyze a wide variety of tumor cells without prior sensitization.

Antioxidants are found to be anticarcinogenic, thereby, protecting the cells against oxidative damage [6]. Bio-normalizer contains antioxidants such as baker's yeast [18,19] and Carica papaya Linn. [19]. Papaya, its major component, was also proven to have antimicrobial and antioxidant properties [20] supplementing its therapeutic action on several cancer growths [21]. Other antioxidative components in Bio-Normalizer include glutathione, hydroxyl-, thiol-, and sulfur-containing amino acids and trace amounts of SOD, catalase and peroxidase [22].

In conclusion, Bio-Normalizer as a nutritional supplement carries the promise of being a potential anti-tumor agent. Its mechanism of action on cancer cells have to be studied in more detail.


References

- Harris, C.C. (1987) Human tissues and cells in carcinogenesis research. Cancer Res., 47, 1-10.
- Santiago, L.A., Osato, J.A., Hiramatsu, M., Edamatsu, R. and Mori, A. (1991) Free radical scavenging action of Bio-catalyzer α.ρ No. 11 (Bio-normalyzer) and its by-product. Free Rad. Biol. Med., 11, 379-383.
- 3. Santiago, L.A., Osato, J.A., Liu, J.and Mori, A. (1993) Age-related increases in superoxide dismutase activity and thiobarbituric acid-reactive substances: Effect of Bio-catalyzer in aged rat brain. Neurochemical Research, 18, 711-717.
- Santiago, L.A., Uno, K., Kishida, T., Miyagawa, F., Osato, J.A. and Mori,
 A. (1994) Effects of Bio-normalizer on serum components and immunological functions in humans. Neurosciences, 20s, 149-152.
- Bernas, G.C., Remo, G.M., Abrigo, R.B., San Miguel, M.M., Osato, J.A., Santiago, L.A., et al. (1993) Antigenotoxic potential of Bio-catalyzer α.p No. 11 (Bio-normalizer) against somatic cell genotoxic agents. Med. Sci. Res., 21, 107-108.
- Sun, Y. (1990) Free radicals, antioxidant enzymes and carcinogenesis.
 Free Rad. Biol. Med., 8, 583-599.
- 7. Flaherty, J.T. (1991) Myocardial injury mediated by oxygen radicals. The Am. J. Med., 91(3C), 79-85.
- 8. Mori, A., Hiramatsu, M., Yokoi, I. and Edamatsu, R. (1990) Biochemical pathogenesis of post-traumatic epilepsy. Pav. J. Biol. Sci. 25, 54-62.
- 9. Breimer, L.H. (1988) Ionizing radiation-induced mutagenesis. Br. J. Cancer, 57, 6 18.
- Schraufstatler, I.U., Hinshaw, D.B., Hyslop, P.A., Spragg, R.G. and Cochrane, C.G. (1986) Oxidant injury of cells: DNA strand-breaks activate polyadenosine diphosphate-ribose polymerase and lead to depletion of nicotinamide adenine dinucleotide. J. Clin. Invest., 77, 1312.
- 11. Cochraine, C.G. (1991) Cellular injury by oxidants. The Am.J. Med., 91(3C), 23-30.
- 12. Clemens, M.R. (1991) Free radicals in chemical carcinogenesis. Klin-Wochenschr, 69, 1123-1134.
- 13. Mossman, T. (1983) Rapid colorimetric assay for cellular growth and

- survival: Application to proliferation and cytotoxic assays. J. Immunol. Methods, 65, 55-63.
- 14. Hansen, M.B., Nielsen, S.E. and Berg K. (1989) Re-examination and further development of a precise and rapid dye method for measuring cell growth/cell kill. J. Immunol. Methods, 119, 203-210.
- Freshney, R.I. (1986) Animal Cell Culture: A Practical Approach, pp. 213-214. Allan R. Liss, Inc., New York.
- Okuda, H., Matsura, Y., Zhou, A., Osato, J.A. and Santiago, L.A. (1993)
 Inhibitory action of Bio-Normalizer on toxohormone-L-induced lipolysis.
 The Clinical Report, 27, 47-56.
- Zhou, A., Matsura, Y. and Okuda, H. (1993) Effect of Bio-catalyzer (Bio-Normalizer) on natural killer cell activity in mice tumor. Clinical Report, 27, 47-48.
- 18. Santiago, L.A., Hiramatsu, M. and Mori, A. (1991) Scavenging action on free radicals and inhibition of lipid peroxidation by baker's yeast. Med. Sci. Res., 19, 867-868.
- 19 Santiago, L.A., Osato, J.A. and Mori, A. (1992) Free radical scavenging action of *Carica papaya* Linn. Neurosciences, 18, 189-194.
- 20. Osato, J.A., Santiago, L.A., Remo, G.M., Cuadra, M.S. and Mori, A. (1993) Antimicrobial and antioxidant activities of unripe papaya. Life Sciences, 53, 1383-1389.
- Quisumbing, E. (1978) Medicinal Plants of the Philippines, pp.632-635.
 Katha Publishing Co., Philippines.
- Osato, J.A., Cuadra, M.S., Santiago, L.A., Mori, A. and Horitsu, H. (1995) Antioxidant defenses of Bio-normalizer. Magnetic Resonance in Medicine, 6, 306-308.

Figure Legend:

- Figure 1 Effects of PBS, Bio-Normalizer (BN) and vincristine sulfate (VS, 20 μg/ml) on HeLa (■), HEp-2 (□), A 549 (□) and SLNI-52 (□). Cytotoxicity of BN and VS on A 549 and SLNI-52 has no significant difference with one another (p > 0.05). BN's IC₅₀ value for each cell line was used.
- Figure 2 Percentage cell viability of HeLa (\Box , p < 0.05), HEp-2 (\bullet , p < 0.05), A 549 (\Box , p < 0.05) and SLNI-52 (\bullet , p < 0.05) cell lines after 4-day exposure to Bio-Normalizer (BN). The points represent the means for 4 determinations.



Test Compounds

Figure 1




Figure 2